



**2002**  
OSISOFT USERS CONFERENCE

**EXPANDING  
THE POWER OF PI**

**MONTEREY CALIFORNIA**



**OSIsoft™**

**Regulatory Compliant *PI* Batch Reporting**

**Designing for the PI Application Framework**

# What is the need ?


■ In 1997 FDA came out with Regulatory Requirement 21 CFR Part 11 – Electronic Records and Signatures.

■ With historians so far the emphasis has been on the ability to store Electronic Records in a secure and auditable manner.

■ But the data is now in a secure 'locked box'. So how can you create batch reports in a compliant manner ?



OSI is now focused on using PI's power to provide compliant reports out of the 'locked box' !







# PI Tools

vs.

# PI Applications

- Broad functionality
- Cross-industry
- No versioning, audit trails
- Unstructured configuration
- User-maintained applications
- Intended for engineering use

- Specific functionality
- Industry-focused
- Integrated security
- Strict configuration control
- Strict version control
- Intended for production use



# Build it and They will Come !!

- It needs to be a fully configurable *PI* batch reporting product.
- Must fully integrate current *PI* automation functionality
  - use the PI MDb, PI API, PI SDK and PBook automation
- It must produce *PI* batch reports without programming:
  - Production Summaries, Exception Reports .....
  - Detailed Batch Activity, Equipment Preparation .....
- Initially targeted at FDA 21 CFR Part 11 reporting compliance using *PI*
- Version 1 has been operational with Biogen since Jan 01.
- Share our plans and design for the *PI* Application Framework


## Customer Batch Report for 750 L Bioreactor

1st Establish

### Batch Summary

Batch Number: PI\_TEST20A  
Product: Test Product  
Batch Started: Tuesday, March 28, 2000 1:45 PM  
Batch Ended: Tuesday, March 28, 2000 5:35 PM  
Batch Duration: 3 hours 50 minutes

User Requirements



### Report Summary

Report Print Date: Thursday, April 19, 2000 4:30 PM  
Report Name: Batch Excursion Report for Bioreactors  
Report Version: Version 1.0  
Report Author: A JFenn  
Report Description: Used to flag any exceptions from normal limits  
Report Template: ReportServer/C:/Reports/Simple Batch Report.rpt  
Report Printed:

### Excursion Summary

750 LBioreactor 15-TK-20A

<u>Tag Name</u>	<u>Limits</u>	<u>Phase</u>	
Vessel Pressure ( psig)	Lower:1 Upper:10	Heating	
	Low: March 28, 2000 1:47 PM		0.0
Vessel Temperature ( deg.C)	Lower:34.5 Upper:36.5	Heating	
	High: March 28, 2000 1:48 PM		37.2
Agitator Speed (RPM)	Lower:20 Upper:40	Heating	
	Low: March 28, 2000 1:57 PM		0.0
{etc.}			



## Process Step Excursion Details

Step Name	Step Start	Step End
Heating 2:10 PM	March 28, 2000 1:45 PM	March 28, 2000

Vessel Pressure (psig)                      Lower Limit : 1                      Upper Limit: 10

Low : March 28, 2000 1:47 PM                      0.0  
Return to Normal: March 28, 2000 1:55 PM

**{low limit violation, with return to normal action}**

Vessel Temperature (deg. C)                      Lower Limit : 34.5                      Upper Limit: 36.5

High : March 28, 2000 1:48 PM                      37.2  
Dissolved Oxygen (%pO2):                      34.3


**{hi limit violation, with log other tag action}**

Agitator Speed (RPM)                      Lower Limit : 20                      Upper Limit: 40



Low : March 28, 2000 1:57 PM                      0.0

Trend: Agitator Speed (RPM)


Agitator Speed



# So how do you currently do this with PI 3.3 ?



Batch Report Viewer



# Batch Report Editor - report building - versioning - tag/limits/actions - release

**InfoBatch - Configurator**

File Edit View Options Reference Windows About

Report Group Report Rep

**Report Explorer**

- DEFAULT REPORT GRO
  - 750 L BioReactor 15-T
  - 2000L BioReactor 15-T
  - 2000L BioReactor 15-T
  - 100L BioReactor 15-T
  - 100L BioReactor 15-T
  - 750 L BioReactor 15-T
  - Version - 1
    - BATCHNAME
      - 15-TK-20
        - Vessel Pre
        - Vessel Tei
        - TREN
        - FIXED
        - T4
        - Agitator Sp
        - TREN
        - FIXED
        - T4
        - Air Overlay
        - Dissolved
        - TREN
        - FIXED
        - T4
        - pH 1
        - TREN
        - FIXED
  - 100L FTS 05-FL-033 C
  - 100L FTS 05-FL-033A
  - Product Hold Tank 40
  - Product Hold Tank 40
  - 750L BioReactor 15-T
  - 750L BioReactor 15-T
  - 2000L BioReactor 15-T

**Report Specifications**

- 15-TK-20
  - Vessel Pressure
  - Vessel Temperature
  - TREND
  - FIXED(HIGH - 37.0: LOW - 35.0)
  - TABLE
  - Agitator Speed
  - TREND
  - FIXED(HIGH - 40: LOW - 20)
  - TABLE
  - Air Overlay
  - Dissolved Oxygen
  - TREND
  - FIXED(HIGH - 120: LOW - 10)
  - TABLE

**Report Configuration**

Report Promote Report State

DEFAULT REPORT GROUP

BATCHREPORT.XSL

18/2000 11:23:27 AM

Cancel

Screen Mode VIEW

1 object(s) selected 79.0KB

Report Explorer

Report Specifications



**InfoBatch - Configurator**

File Edit View Options

Report Group Rep

**Item Commands**

Copy Item Create Item Edit Item Archive Item

**Create/Edit Item**

DataSource: PI

Source: BPU108:FIC-1149A/PID1/PV.CV Version: 1

Item Name: AVE FLOW RATE

Display Precision: 0.0 Uom: State: COMMITTED

Commit Cancel

**Additional Screens**

Limit Maintenance Action Maintenance Maintenance

**Item-Limit Settings**

DataSource: Limit Type: Limit Class: High: Target: Low:

Create Date: Created By:

Hysteresis:

Associated Actions:

**Item-Action Settings**

Action Name: CALC Action Type: REGULAR Trend Type: Calculation: ARCAVERAGE

Offset: 30 Offset Operator: Increase | NoOffset

Interval:

RCOVERAGE  
Increase | NoOffset


Screen Mode VIEW

UserName: admin

1 object(s) selected

# Highlights of Report Editor

- **Supports PI Module Db**
- Multi level reports
  - Batch, unitop, operation procedure, phase etc...
  - Additional hierarchies
- Copy/Cut/Paste
- Automatic versioning and auditing
- Intelligent pop up menus
- Intelligent use of PI data
- Tag string parsing
- Tag groupings
- Detailed administration functions
  - Multiple Limit types
  - Define product types
  - Registry settings



# Batch Report Generator - batch . unit . product search - report generation - report print

InfoBatch - [15-TK-80B Batch Report - Version 1]

File Edit View

Batches Units

15-TK-20  
15-TK-20A

**BIOGEN**

**Batch Number : CC2-01-002**

---

**Report Summary**

Report Generation Date: **23-Mar-2001 11:41:33 AM** Report Name: **15-TK-80B Batch Report**  
Report Version: **1** Report Author: **T. Fenn**  
Report Printed:  
Report Description: **15-TK-80B SUMMARY REPORT**  
Report Template: **BATCHREPORT.XSL**

---

**Unit Summary**

Unit	Start Time	End Time	Duration
15-TK-80B	19-Mar-2001 11:07:19 AM	IN PROGRESS	12 Hrs 52 Mins

CIP:05-FL-001A

Vessel Pressure Trend


02/19/2001 1:26:44 PM : 03/23/01 2:43:17 AM

Batch Explorer

Browser Report Viewer



# Exception Report showing embedded trends



# Highlights of Report Generator




- **PI Module DB support of batch searching**
- Batch Tree View
- Generate Reports on any SP88 level
- Protected trends
- Multiple tags on trends
- Additional security of embedded Explorer
- Printout only allowed
- Auto pagination of report printouts
- XSL Style sheets embedded in executable
- Deletion of temporary files


# Why Application Framework?

- Infrastructure, Infrastructure, Infrastructure
  - Elements (Equipment)
  - Operating Environment
  - Extensible
  - Components
  - Tighter Application Integration






# PI Application Framework Server




Batch Report Editor


## PI AF PLUG INS




XML




XSL Style Sheets




Batch Report Viewer



# PI Application Framework



# Report Generation Server





# Element/Models

- Equipment Extensions
  - Specific Equipment Usage
  - Product Type
- Report Structure
  - Report Summary Information
  - Report Section Information
  - Specific Limits and Actions


# Plug-Ins

- Batch Time Rule
  - Provide Batch Functionality to Plug-Ins
- Core AF Plug-Ins w/ VBJ Plug-Ins
  - Perform Specific Service, I.e. CheckLimits, Associated Tag Information
- Report Generation Plug-Ins
  - Report Generation Infrastructure
  - Interpret report structure
  - Provide Context information to Core AF Plug-ins

# VBJ Plug-Ins

- Process Event Analysis Rules
  - Limit Check
  - Ramping Check
  - State Change Check
  - Event Monitoring
  - Material Charge Information

# BATCH REPORTING ARCHITECTURE USING THE PI-AF and VBJ Plug-Ins





# Benefits And Future Functionality

- Types of Report
  - Equipment Based Reporting
  - Procedure Based Reporting
  - Standard Reports (Shift, Daily, etc.)
- Future Functionality
  - Report Output/Export
  - Report Approval Process
  - Report Verification
  - Integration w/ ICE