



# PI System in a Lab-Based Course: Bio Engineering at MIT

Presented by **Erica Trump Alexander Kendrick** 






# **Chemical - Biological Engineering Laboratory**

 Revamped, project-based version of Chemical Engineering Unit Operations

- Course objectives include:
  - Hands-on experience with bioprocesses
  - Experimental design
  - Exposure to industry-standard tools



## The Life Sciences Story at OSIsoft

















GILEAD



















#### More reasons for PI


- Real-time monitoring and analysis
- Less data pushing, more time spent on scientific insights and design
- Students pool results
- Retain knowledge from past semesters
- Remote monitoring ability
- Integration with tools such as MATLAB

#### Classroom rollout

- System setup and education of lab staff
- Start small in pilot semester
- Introduce students to PI Visualization Tools
  - Workbooks with step-by-step instructions, tie-ins to OSIsoft Learning YouTube channel, and student exercises
  - Training databases with simulated Unit Operations data

#### Classroom rollout

- Assignment 1
  - PI DataLink to analyze simulated process data
- Assignment 2
  - PI ProcessBook to prepare bioreactor displays to be used throughout course



An example student display

#### Classroom rollout

- Assignment 3
  - Analyze real-time results from experiments
  - Students pool results
    - Explore effects of a wide range of independent variables



Student interaction during a mass transfer experiment



- Outline:
- Improve the Students Understanding Bioreactors with PI
- New Lesson to the Syllabus
  - Integration of Outside Data Analysis

# Working PI into the Syllabus: Future Plans

## Improving the Students Understanding of Bioreactors with PI

Use previous PI data from early semesters to educate the students to use PI

#### Benefits:

- Project relevant teaching
- Increased student motivation
- Improved experimental design



## **New lesson to the Syllabus**

- The new lesson was tested with incoming researchers in the Hamel laboratory in the Spring term (IAP).
- The goal was to not only have the students present data as they became accustomed to PI but to also learn to interpret data.
- In this, we created case studies to get students to begin to think about their reactors and make the connections of what they data is communicating
  - Event frames were implemented to bookmark experimental runs to make finding experiments easier



# **Example of the New Teaching Style**

Students are given an event and the desired parameters that the reactor should keep throughout the process as well as reactor type.

Example parameters:

Reactor type: Batch


Reactor working volume: 1.5L

pH= 7 T=30C

Stir rate: 300rpm


Dissolved Oxygen: not

controlled



Students will be encouraged to make their own displays and access odd occurrences in the trends.

The students will notice pH is amiss (pH≠ 7) and will investigate acid and base addition values



PI Coresight H New Display 4 NewDisplay9 Undo Acid and Base Trend Values > working volume Home ▶ HamelLab ▶ BioRtr10 Q BioRtr10|OD | O BioRtr10|Dissolved Oxyger BioRtr10|pH 0.055658 49.091 2,309 7.0045 pH Exhaust Oxygen 0.064 90 500 2,500 12 OXC RMS flow RQ. Process Values 0.06 Acid Pump Total AntiFoam Pump T ... 0.058 Base Pump Total Dissolved Oxygen 0.056 PH pH Stir Rate 0.054 Substrate Pump Tr Temperature 0.052 System Communication U 0.05 10 0 Lifebeat 5d 11/16 3:41 PM +2d 23h 11/13/2015 3:45:00 PM 11/18/2015 3:45:00 PM 11/13/2015 3:45:00 PM 1h 8h 1d 1w 1mo & Now 11/18/2015 3:45:00 PM Search Related Assets/Events (12/

Oct 4 , , , Oct 11 , , Oct 18 , , , Oct 25 , , , Noy 1 , , , Noy 8 , , , , , Noy 15 , , ,

Students will see that the acid and base additions are illogical and diagnose the problem.

# **Integrating Outside Data Analysis Tools**


# We plan to introduce MatLab data manipulation to the current syllabus

- Students are already familiar with Matlab by this point in their studies
- A program, the PI-to-MATLAB Utility, supplied by OSIsoft already exists and is easily accessible
- The additional computing capacity will allow the students complete the carbon mass balance to analysis cellular metabolism

## The PI-to-MATLAB Utility



The Pi-To-MatLab shortcut directly opens MatLab and allows you to access the PI server is credentials are acceptable



#### **Contact Information**

Erica Trump, PhD

etrump@osisoft.com
Instructional Systems Designer
OSIsoft

Alexander Kendrick
<a href="mailto:ankendri@mit.edu">ankendri@mit.edu</a>
Laboratory Assistant/Researcher
MIT

Jean-François Hamel, PhD jhamel@mit.edu Principal Investigator MIT







감사합니다

Danke

谢谢

Gracias

Thank You

ありがとう

Спасибо

Obrigado



Merci